

## Product datasheet

# Anti-Netrin 4 antibody [191633] ab50268

### 製品の概要

製品名	Anti-Netrin 4 antibody [191633]
製品の詳細	Mouse monoclonal [191633] to Netrin 4
由来種	Mouse
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	Recombinant full length protein (Human)

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: 5% Trehalose, PBS
精製度	Protein G purified
ポリモノ	モノクローナル
クローン名	191633
ミエローマ	NS0
アイソタイプ	IgG2a

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab50268** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

ELISA

### 追加情報

Direct ELISA: Use at a concentration of 0.5 - 1 µg/ml.

WB: Use at a concentration of 1 - 2 µg/ml. Predicted molecular weight: 70 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

機能	May play an important role in neural, kidney and vascular development. Promotes neurite elongation from olfactory bulb explants.
組織特異性	Expressed in kidney, spleen, mammary gland, aorta, heart, ovary, prostate and fetal spleen.
配列類似性	Contains 3 laminin EGF-like domains. Contains 1 laminin N-terminal domain. Contains 1 NTR domain.
細胞内局在	Secreted > extracellular space > extracellular matrix > basement membrane. Major component.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors