

## Product datasheet

# Anti-Cytokeratin 15 antibody ab51055

### 画像数 1

#### 製品の概要

製品名	Anti-Cytokeratin 15 antibody
製品の詳細	Rabbit polyclonal to Cytokeratin 15
アプリケーション	<b>適用あり:</b> WB, IP <b>適用なし:</b> Flow Cyt or ICC/IF
種交差性	<b>交差種:</b> Human <b>非交差種:</b> Mouse, Rat
免疫原	Synthetic peptide corresponding to residues on the N terminus of human Cytokeratin 15.
ポジティブ・コントロール	WB: A431 cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	PBS 49%, Sodium azide 0.01%, Glycerol 50%, BSA 0.05%
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab51055** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
IP		

追加情報	IP: 1/40. WB: 1/20000. Detects a band of approximately 49 kDa (predicted molecular weight: 49 kDa).  Is unsuitable for Flow Cyt, IHC or ICC/IF.
------	--

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


### 組織特異性

Expressed in a discontinuous manner in the basal cell layer of adult skin epidermis, but continuously in the basal layer of fetal skin epidermis and nail. Also expressed in the outer root sheath above the hair bulb in hair follicle (at protein level). Expressed homogeneously in all cell layers of the esophagus and exocervix, but detected in the basal cell layer only of oral mucosa, skin and in the basal plus the next two layers of the suprabasal epithelium of the palate.

### 配列類似性

Belongs to the intermediate filament family.

## Anti-Cytokeratin 15 antibody 画像



Anti-Cytokeratin 15 antibody (ab51055) at  
1/20000 dilution + A431 cell lysate at 10 µg

### Secondary

Goat anti-Rabbit HRP labeled. at 1/2000  
dilution

**Predicted band size** : 49 kDa

**Observed band size** : 49 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

