

Product datasheet

Anti-NK1.1 antibody [PK136] (PE/Cy7[®]) ab51239

画像数 1

製品の概要

製品名	Anti-NK1.1 antibody [PK136] (PE/Cy7 [®])
製品の詳細	Mouse monoclonal [PK136] to NK1.1 (PE/Cy7 [®])
由来種	Mouse
標識	PE/Cy7 [®] . Ex: 496nm, Em: 774nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	CD mouse spleen and bone marrow cells enriched for NK-1+ cells
特記事項	This product or portions thereof is manufactured under license from Carnegie Mellon University under U.S. Patent Number 5,268,486 and related patents. Cy and CyDye are trademarks of GE Healthcare Limited.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.09% Sodium azide Constituents: 0.14% Monobasic dihydrogen sodium phosphate, 0.87% Sodium chloride, 0.1% Gelatin
	Also contains stabilising agent.
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	PK136
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab51239** in the following tested applications.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors