

Product datasheet

Anti-GANP antibody ab53083

画像数 1

製品の概要

製品名	Anti-GANP antibody
製品の詳細	Rabbit polyclonal to GANP
特異性	GANP antibody detects endogenous levels of total GANP protein.
アプリケーション	適用あり: IHC-P, ELISA
種交差性	交差種: Mouse, Human
免疫原	Synthesized peptide derived from human GANP.
ポジティブ・コントロール	Human breast carcinoma tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab53083** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.
ELISA		1/5000.

ターゲット情報

機能	May be involved in the nuclear localization pathway of MCM3.
----	--


配列類似性

Belongs to the SAC3 family.

細胞内局在

Cytoplasm. Nucleus.

Anti-GANP antibody 画像



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using ab53083 at 1/50 dilution, with and without blocking peptide (ab53083 immunogen)

No blocking peptide

With blocking peptide

Immunohistochemistry (Paraffin-embedded sections)

- GANP antibody (ab53083)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors