

## Product datasheet

# Anti-Dlx2 antibody ab54673

### 画像数 2

#### 製品の概要

製品名	Anti-Dlx2 antibody
製品の詳細	Mouse monoclonal to Dlx2
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Mouse, Human
免疫原	Recombinant fragment, corresponding to amino acids 1-110 of Human Dlx2

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab54673** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
IHC-P		

#### 追加情報

WB: Use at a concentration of 1-5 µg/ml.  
IHC-P: Use at a concentration of 5 µg/ml.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


機能	Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.
配列類似性	Belongs to the distal-less homeobox family. Contains 1 homeobox DNA-binding domain.
細胞内局在	Nucleus.

## 画像



IHC-P - Dlx2 antibody (ab54673)

Dlx2 antibody (ab54673) used in immunohistochemistry at 5ug/ml on formalin fixed and paraffin embedded human stomach.



Western blot - Dlx2 antibody (ab54673)

**Predicted band size : 34 kDa**  
Dlx2 antibody (ab54673) at 1ug/lane +  
NIH/3T3 cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors