

Product datasheet

Anti-Integrin beta 5 antibody ab55359

2 References 画像数 1

製品の概要

製品名	Anti-Integrin beta 5 antibody
製品の詳細	Mouse monoclonal to Integrin beta 5
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment: TASFEVSLEA RSCPSRHTEH VFALRPVGFR DSLEVGVTYN CTGCGSVGLE PNSARCNGSG TYWCGLCCECS PGYLGTRCEC QDGENQSVYQ NLCREA, corresponding to amino acids 421-517 of Human Integrin beta 5 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

特記事項

Abcam is committed to meeting high standards of ethical manufacturing and has decided to discontinue this product by June 2019 as it has been generated by the ascites method. We are sorry for any inconvenience this may cause. We would recommend antibody [ab184312](#) as a replacement.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab55359** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 88 kDa.

## ターゲット情報

機能	Integrin alpha-V/beta-5 is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand.
配列類似性	Belongs to the integrin beta chain family. Contains 1 PSI domain. Contains 1 VWFA domain.
細胞内局在	Membrane.

## 画像



Integrin beta 5 antibody (ab55359) at 1ug/lane  
+ HeLa cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors