

Product datasheet

Anti-Methionyl tRNA synthetase antibody ab55455

画像数 2

製品の概要

製品名	Anti-Methionyl tRNA synthetase antibody
製品の詳細	Mouse monoclonal to Methionyl tRNA synthetase
由来種	Mouse
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	Recombinant fragment: LFQKLENDQI ESLRQRFGGG QAKTSPKPAV VETVTTAKPQ QIQALMDEV T KQGNVRELK AQKADKNEVA AEVAKLLDLK KQLAVAEGKP PEAPKGKKK, corresponding to amino acids 801-900 of Human Methionyl tRNA synthetase Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab55455** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 101 kDa.


アプリケーション	Abreviews	特記事項
----------	-----------	------

ICC/IF		Use a concentration of 10 µg/ml.
--------	--	----------------------------------

ターゲット情報


配列類似性	<p>Belongs to the class-I aminoacyl-tRNA synthetase family.</p> <p>Contains 1 GST C-terminal domain.</p> <p>Contains 1 WHEP-TRS domain.</p>
細胞内局在	Cytoplasm.

画像



Western blot - Methionyl tRNA synthetase antibody (ab55455)

Methionyl tRNA synthetase antibody (ab55455) at 1ug/lane + HepG2 cell lysate at 25ug/lane.



Immunocytochemistry/ Immunofluorescence- Methionyl tRNA synthetase antibody(ab55455)

ab55455 at 10 ug/ml staining Methionyl tRNA synthetase in human Hella cells by Immunocytochemistry/ Immunofluorescence.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors