

Product datasheet

Anti-PNMA1 antibody ab56496

画像数 2

製品の概要

製品名	Anti-PNMA1 antibody
製品の詳細	Mouse monoclonal to PNMA1
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment: KAALLELTGA VDYAAIPREM PGKGGVWKVL FKPPTSDAEF LERLHLFLAR EGWTVQDVAR VLGFNPTPT PGPEMPAEML NYILDNVIQP LVESIWKYR, corresponding to amino acids 62-161 of Human PNMA1 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56496** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
IHC-P		

追加情報

WB: Use at a concentration of 1-5 µg/ml.

IHC-P: Use at a concentration of 3 µg/ml.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

### 組織特異性

Testis and brain specific. In some patients suffering from cancers, it is also specifically expressed by the paraneoplastic tumor cells.


### 配列類似性

Belongs to the PNMA family.

### 細胞内局在


Nucleus > nucleolus. In tumor cells, it is cytoplasmic.

## Anti-PNMA1 antibody 画像



IHC-P - PNMA1 antibody (ab56496)

PNMA1 antibody (ab56496) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human stomach.



Western blot - PNMA1 antibody (ab56496)

**Predicted band size : 40 kDa**

PNMA1 antibody (ab56496) at 1ug/lane + IMR-32 cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors