

Product datasheet

Anti-PNMA1 antibody ab56496

画像数 2

製品の概要

製品名	Anti-PNMA1 antibody
製品の詳細	Mouse monoclonal to PNMA1
由来種	Mouse
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment: KAALLELTGA VDYAAIPREM PGKGGVWKVL FKPPTSDAEF LERLHLFLAR EGWTVQDVAR VLGFNPTPT PGPEMPAEML NYILDNVIQP LVESWYKR, corresponding to amino acids 62-161 of Human PNMA1 Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56496** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 40 kDa.


アプリケーション	Abreviews	特記事項
----------	-----------	------

IHC-P		Use a concentration of 3 µg/ml.
-------	--	---------------------------------

ターゲット情報


組織特異性	Testis and brain specific. In some patients suffering from cancers, it is also specifically expressed by the paraneoplastic tumor cells.
配列類似性	Belongs to the PNMA family.
細胞内局在	Nucleus > nucleolus. In tumor cells, it is cytoplasmic.

画像



PNMA1 antibody (ab56496) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human stomach.

IHC-P - PNMA1 antibody (ab56496)



PNMA1 antibody (ab56496) at 1ug/lane + IMR-32 cell lysate at 25ug/lane.

Western blot - PNMA1 antibody (ab56496)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors