

Product datasheet

Anti-PNMA1 antibody ab56496

画像数 2

製品の概要

製品名	Anti-PNMA1 antibody
製品の詳細	Mouse monoclonal to PNMA1
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment: KAALLELTGA VDYAAIPREM PGKGGVWKVL FKPPTSDAEF LERLHLFLAR EGWTVQDVAR VLGFNQPTPT PGPEMPAEML NYILDNVIQP LVESWYKR, corresponding to amino acids 62-161 of Human PNMA1 Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56496** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
IHC-P		

追加情報

WB: Use at a concentration of 1-5 µg/ml.

IHC-P: Use at a concentration of 3 µg/ml.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

組織特異性

Testis and brain specific. In some patients suffering from cancers, it is also specifically expressed by the paraneoplastic tumor cells.


配列類似性

Belongs to the PNMA family.

細胞内局在


Nucleus > nucleolus. In tumor cells, it is cytoplasmic.

画像



IHC-P - PNMA1 antibody (ab56496)

PNMA1 antibody (ab56496) used in immunohistochemistry at 3µg/ml on formalin fixed and paraffin embedded human stomach.



Western blot - PNMA1 antibody (ab56496)

Predicted band size : 40 kDa

PNMA1 antibody (ab56496) at 1 µg/lane + IMR-32 cell lysate at 25 µg/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors