

Product datasheet

Anti-Syntaxin 8 antibody ab56667

画像数 1

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Syntaxin 8 antibody |
| 製品の詳細 | Mouse monoclonal to Syntaxin 8 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant full length protein, corresponding to amino acids 1-237 of Human Syntaxin 8 |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Preservative: None PBS, pH 7.2 |
| 精製度 | Protein G purified |
| ポリ/モノ | モノクローナル |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56667** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 27 kDa. |

ターゲット情報

| | |
|--------------|---|
| 機能 | Vesicle trafficking protein that functions in the early secretory pathway, possibly by mediating retrograde transport from cis-Golgi membranes to the ER. |
| 組織特異性 | Highly expressed in heart. Also found in brain, kidney, liver, lung, placenta, skeletal muscle, spleen and pancreas. |
| 配列類似性 | Belongs to the syntaxin family. Contains 1 t-SNARE coiled-coil homology domain. |
| 翻訳後修飾 | Ubiquitinated by HECTD3. |
| 細胞内局在 | Membrane. Preferentially associated with the early endosome. To lesser extends, also present in late endosome, the plasma membrane and coated pits. |

画像



Syntaxin 8 antibody (ab56667) at 1ug/lane +
A-431 cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors