

## Product datasheet

# Anti-SPAG6 antibody ab56778

★★★★★ 1 Abreviews [画像数 2](#)

### 製品の概要

製品名	Anti-SPAG6 antibody
製品の詳細	Mouse monoclonal to SPAG6
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Human
免疫原	Recombinant fragment: MSQRQVLQVF EQYQKARTQF VQMVAELATR PQNIETLQNA GVMSLLRLL LDVVPTIQQT AALALGRLAN YNDDLAEAVV KCDILPQLVY SLAEQNRFYK , corresponding to amino acids 1-101 of Human SPAG6 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG3
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56778** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★★★★★	Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 55 kDa.

## ターゲット情報


---

機能	Important for structural integrity of the central apparatus in the sperm tail and for flagellar motility.
組織特異性	Highly expressed in testis.
配列類似性	Contains 8 ARM repeats.
細胞内局在	Cytoplasm > cytoskeleton. Cell projection > cilium > flagellum. Associated with microtubules. Detected on the sperm flagellum.

---


## 画像

---



Western blot - SPAG6 antibody (ab56778)

SPAG6 antibody (ab56778) at 1ug/lane + A-431 cell lysate at 25ug/lane.



Western blot - SPAG6 antibody (ab56778)  
This image is courtesy of an anonymous abreview.

Anti-SPAG6 antibody (ab56778) at 1/500 dilution + whole tissue lysate prepared from mouse lung at 80 µg

**Secondary**

Sheep anti-mouse conjugated to HRP at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 55 kDa

**Observed band size:** 55 kDa

**Exposure time:** 20 minutes

Primary antibody incubated for 2 hours at 25°C

10% PAGE. Denaturing conditions

Blocked with 1% milk for 16 hours at 4°C.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors