

## Product datasheet

# Anti-PTPLAD1 antibody ab57143

1 References [画像数 1](#)

### 製品の概要

製品名	Anti-PTPLAD1 antibody
製品の詳細	Mouse monoclonal to PTPLAD1
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to amino acids 1-114 of Human PTPLAD1

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG1
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab57143** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB		Use a concentration of 1 - 5 µg/ml. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.
----	--	---

## ターゲット情報


---

**関連性** PTPLAD1 is a component of Rac1-signaling pathways leading to the modulation of gene expression. It potentiates different effects of the small GTPase Rac1, such as c-Jun N-terminal kinase activation and transcriptional activity of nuclear factor kappaB (NF-kappaB).It plays a crucial role in HCV RNA replication and the propagation of JFH1 virus through interaction with viral and host proteins.

**細胞内局在** Membrane; Multi-pass membrane protein

## 画像

---



Western blot against tagged recombinant protein immunogen using ab57143 PTPLAD1 antibody at 1ug/ml. Predicted band size of immunogen is 39 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors