

Product datasheet

Anti-MFRP antibody ab57797

画像数 1

製品の概要

製品名	Anti-MFRP antibody
製品の詳細	Mouse monoclonal to MFRP
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment: NTTAFPNMV GMTQEEVVE VLSGYKSLTS LPCYQHFRRL LCGLLVPRCT PLGSVLPPCR SVCQEAHQ QSGLALLGTP WPFNCNRLPE AADLEACAQP , corresponding to amino acids 480-580 of Human MFRP Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2b
軽鎖の種類	lambda

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab57797** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

追加情報

WB: Use at a concentration of 1-5 µg/ml.

This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報


関連性

Membrane frizzled-related protein, MFRP, may have a role in eye development and may be necessary for photoreceptor maintenance. Mutations of membrane-type frizzled-related protein are associated with inherited retinal degenerations including syndromic forms of microphthalmos/retinal degeneration, a rare autosomal recessive disorder of eye development characterized by extreme hyperopia and small functional eyes, posterior microphthalmia with retinitis pigmentosa, foveoschisis, and optic disk drusen.

細胞内局在

Cell Membrane; single-pass type II membrane protein

画像



Western blot against tagged recombinant protein immunogen using ab57797 MFRP antibody at 1ug/ml. Predicted band size of immunogen is 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors