

Product datasheet

Anti-TAL2 antibody ab58236

画像数 1

製品の概要

製品名	Anti-TAL2 antibody
製品の詳細	Mouse monoclonal to TAL2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment: IPTHPPDKKL SKNETLRLAM RYNFLVKVL GEQSLQQTGV AAQGNILGLF PQGPHLPGLE DRTLLENYQV PSPGPSHHIP , corresponding to amino acids 29-109 of Human TAL2 Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab58236** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

アプリケーション	Abreviews	特記事項
----------	-----------	------


WB Use a concentration of 1 - 5 µg/ml.
 This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

ターゲット情報

関連疾患 Note=A chromosomal aberration involving TAL2 may be a cause of some T-cell acute lymphoblastic leukemia (T-ALL). Translocation t(7;9)(q34;q32) with TCRB.

配列類似性 Contains 1 basic helix-loop-helix (bHLH) domain.

画像



Western blot against tagged recombinant protein immunogen using ab58236 TAL2 antibody at 1ug/ml. Predicted band size of immunogen is 35 kDa

Western blot - TAL2 antibody (ab58236)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors