

Product datasheet

Anti-LRRN5 antibody ab58276

画像数 1

製品の概要

製品名	Anti-LRRN5 antibody
製品の詳細	Mouse monoclonal to LRRN5
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment: RVQETHPYHI LLSWVTPPNT VSTNLTWSSA SSLRGQGATA LARLPRGTHS YNITRLQAT EYWACLQVAF ADAHTQLACV WARTKEATSC HRALGDRPG, corresponding to amino acids 532-631 of Human LRRN5 Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab58276** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

追加情報

WB: Use at a concentration of 1-5 µg/ml.

This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性 LRRN5 (Leucine rich repeat neuronal 5) is a transmembrane protein belonging to the leucine-rich repeat (LRR) superfamily. It has homology with other proteins that function as cell-adhesion molecules or as signal transduction receptors. It is overexpressed in malignant gliomas.

細胞内局在 Membrane; Single-pass type I membrane protein.

画像



Western blot against tagged recombinant protein immunogen using ab58276 LRRN5 antibody at 1ug/ml. Predicted band size of immunogen is 37 kDa

Western blot - LRRN5 antibody (ab58276)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors