

Product datasheet

Anti-OAS2 antibody ab60091

画像数 1

製品の概要

製品名	Anti-OAS2 antibody
製品の詳細	Rabbit polyclonal to OAS2
由来種	Rabbit
アプリケーション	適用あり: ELISA, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide (10-30 aa in length) at the C-term of last 50 aa of human OAS2 conjugated to KLH
ポジティブ・コントロール	Human breast carcinoma tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Protein G purified
特記事項 (精製)	ab60091 is purified through a protein G column and eluted out with both high and low pH buffers and neutralized immediately after elution then followed by dialysis against PBS.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab60091** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ELISA		1/1000.

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/100.

ターゲット情報

機能	May play a role in mediating resistance to virus infection, control of cell growth, differentiation, and apoptosis.
配列類似性	Belongs to the 2-5A synthase family.
細胞内局在	Mitochondrion. Nucleus. Microsome. Endoplasmic reticulum. Associated with different subcellular fractions such as mitochondrial, nuclear, and rough/smooth microsomal fractions.

画像



ab60091 at 1/50 - 1/100 dilution, staining OAS2 in human breast carcinoma by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue, followed by peroxidase-conjugated secondary antibody and AEC staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - OAS2 antibody (ab60091)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors