


Product datasheet

Anti-DAP Kinase 2 (phospho S318) antibody ab61042

画像数 1

製品の概要

製品名	Anti-DAP Kinase 2 (phospho S318) antibody
製品の詳細	Rabbit polyclonal to DAP Kinase 2 (phospho S318)
特異性	Detects endogenous levels of DAP Kinase 2 only when phosphorylated at serine 318.
アプリケーション	適用あり: ELISA, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse 
免疫原	Synthetic phosphopeptide derived from human DAP Kinase 2 around the phosphorylation site of serine 318 (K-L-S ^P -F-S).
ポジティブ・コントロール	Skeletal muscle

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
精製度	Immunogen affinity purified
特記事項 (精製)	Affinity purified from rabbit antiserum by affinity chromatography using epitope specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatography using non-phosphopeptide corresponding to the phosphorylation site.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab61042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/5000.
IHC-P		1/50 - 1/100.

ターゲット情報

機能	Calcium/calmodulin-dependent serine/threonine kinase which acts as a positive regulator of apoptosis.
組織特異性	Ubiquitously expressed in all tissue types examined. High levels in heart, lung and skeletal muscle.
配列類似性	Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. DAP kinase subfamily. Contains 1 protein kinase domain.
細胞内局在	Cytoplasm.

画像



1/50-1/100 ab61042 staining phosphorylated DAP Kinase 2 in human skeletal muscle tissue, in the absence or presence of the immunizing phosphopeptide.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors