


Product datasheet

Anti-GDNF Receptor alpha 1 antibody ab62123

画像数 1

製品の概要

製品名	Anti-GDNF Receptor alpha 1 antibody
製品の詳細	Rabbit polyclonal to GDNF Receptor alpha 1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human 交差が予測される動物種: Cow 
免疫原	Synthetic peptide: TALRVKKNKPLGPAG , corresponding to internal sequence amino acids 369-382 of Human GDNF Receptor alpha 1 Run BLAST with Run BLAST with
ポジティブ・コントロール	Membrane fraction lysates of human brain, liver, kidney and spleen.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Protein A purified
特記事項(精製)	Purified IgG
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab62123** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500. Detects a band of approximately 51 kDa (predicted molecular weight: 51 kDa).

ターゲット情報

機能	Receptor for GDNF. Mediates the GDNF-induced autophosphorylation and activation of the RET receptor.
配列類似性	Belongs to the GDNFR family.
細胞内局在	Cell membrane.

画像



All lanes : Anti-GDNF Receptor alpha 1 antibody (ab62123) at 1/500 dilution

Lane 1 : Crude membrane fraction of human brain

Lane 2 : Crude membrane fraction of human liver

Lane 3 : Crude membrane fraction of human kidney

Lane 4 : Crude membrane fraction of human spleen

Lysates/proteins at 15 µg per lane.

Secondary

All lanes : Anti-rabbit IgG HRP

Predicted band size: 51 kDa

Observed band size: 51 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors