

Product datasheet

Anti-CNOT4 antibody ab63028

画像数 1

製品の概要

製品名	Anti-CNOT4 antibody
製品の詳細	Rabbit polyclonal to CNOT4
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	Synthetic peptide from the C-terminal region of human CNOT4 Isoform A, conjugated to KLH.
ポジティブ・コントロール	CEM cell lysates.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Protein G purified
特記事項 (精製)	ab63028 was purified through a protein G column, eluted with high and low pH buffers and neutralized immediately after elution, followed by dialysis against PBS.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab63028** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/50 - 1/100. Detects a band of approximately 60 kDa (predicted molecular weight: 64 kDa).
ELISA		1/1000.

ターゲット情報

機能	Has E3 ubiquitin ligase activity. The CCR4-NOT complex functions as general transcription regulation complex.
パスウェイ	Protein modification; protein ubiquitination.
配列類似性	Contains 1 C3H1-type zinc finger. Contains 1 RING-type zinc finger. Contains 1 RRM (RNA recognition motif) domain.
翻訳後修飾	Autoubiquitinated. Phosphorylated upon DNA damage, probably by ATM or ATR.
細胞内局在	Cytoplasm. Nucleus.

画像



Anti-CNOT4 antibody (ab63028) at 1/60 dilution + CEM cell line lysates at 35 µg/ml

Predicted band size: 64 kDa

Observed band size: 60 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors