

Product datasheet

Anti-Ninjurin2 antibody ab67140

画像数 1

製品の概要

製品名	Anti-Ninjurin2 antibody
製品の詳細	Mouse polyclonal to Ninjurin2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length Ninjurin2 protein (Human)
ポジティブ・コントロール	Ninjurin2-transfected 293T cells lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67140** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Detects a band of approximately 16 kDa (predicted molecular weight: 16 kDa).

ターゲット情報

機能	Homophilic cell adhesion molecule that promotes axonal growth. May play a role in nerve regeneration and in the formation and function of other tissues.
組織特異性	Widely expressed. In adult, higher expression in the bone marrow and peripheral blood lymphocytes, medium in the lung, lymph node, thyroid, uterus, thymus, spleen, prostate and skeletal muscle, lower in the liver, placenta, brain, heart and kidney. In embryo, higher expression in the thymus, heart and liver, lower in the spleen, lung, brain and kidney.
配列類似性	Belongs to the ninjurin family.
細胞内局在	Membrane.

画像



All lanes : Anti-Ninjurin2 antibody (ab67140)
at 1/500 dilution

Lane 1 : Ninjurin2-transfected 293T cells
lysate

Lane 2 : non-transfected 293T cells lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat anti-mouse IgG (H&L) HRP at
1/2500 dilution

Predicted band size: 16 kDa

Observed band size: 16 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors