

Product datasheet

Anti-C8orf72 antibody ab67209

画像数 1

製品の概要

製品名	Anti-C8orf72 antibody
製品の詳細	Mouse polyclonal to C8orf72
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein (Human)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67209** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 51 kDa (predicted molecular weight: 41 kDa).

ターゲット情報

機能	May be involved in tumor progression.
----	---------------------------------------

組織特異性


Detected in thyroid, spleen and testis, and at lower levels in stomach, spinal cord, lymph node, trachea, adrenal gland, prostate, ovary and intestine.

配列類似性

Belongs to the FAM110 family.

細胞内局在

Cytoplasm. Cytoplasm > cytoskeleton > centrosome.

画像

All lanes : Anti-C8orf72 antibody (ab67209)
at 1/500 dilution

Lane 1 : C8orf72 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Predicted band size: 41 kDa

Observed band size: 51 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors