

## Product datasheet

# Anti-ZNF434 antibody ab67291

### 画像数 1

#### 製品の概要

|              |                                     |
|--------------|-------------------------------------|
| 製品名          | Anti-ZNF434 antibody                |
| 製品の詳細        | Mouse polyclonal to ZNF434          |
| 由来種          | Mouse                               |
| アプリケーション     | 適用あり: WB                            |
| 種交差性         | 交差種: Human                          |
| 免疫原          | Full length ZNF434 protein (Human)  |
| ポジティブ・コントロール | ZNF434 transfected 293T cell lysate |

#### 製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Preservative: None<br>Constituents: Whole serum   |
| 精製度    | Whole antiserum   |
| ポリ/モノ  | ポリクローナル   |
| アイソタイプ | IgG   |

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67291** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| WB       |           | 1/500 - 1/1000. Detects a band of approximately 28 kDa (predicted molecular weight: 55 kDa). |

#### ターゲット情報

|              |  |
|--------------|--|
| <b>機能</b>    | May be involved in transcriptional regulation.   |
| <b>配列類似性</b> | Belongs to the krueppel C2H2-type zinc-finger protein family.<br>Contains 6 C2H2-type zinc fingers.<br>Contains 1 KRAB domain. |
| <b>細胞内局在</b> | Nucleus.   |

## 画像



Western blot - ZNF434 antibody (ab67291)

**All lanes :** Anti-ZNF434 antibody (ab67291)  
at 1/500 dilution

**Lane 1 :** ZNF434 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg/ml per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP  
Conjugate at 1/2500 dilution

**Predicted band size:** 55 kDa

**Observed band size:** 28 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors