

## Product datasheet

# Anti-NPM2 antibody ab67403

★★★★☆ 1 Abreviews 1 References 画像数 1

### 製品の概要

製品名	Anti-NPM2 antibody
製品の詳細	Mouse polyclonal to NPM2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-214 of Human NPM2
ポジティブ・コントロール	Transfected 293T cells

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67403** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 34 kDa (predicted molecular weight: 24 kDa).

### ターゲット情報

**機能** Core histones chaperone involved in chromatin reprogramming, specially during fertilization and

early embryonic development. Probably involved in sperm DNA decondensation during fertilization.

**配列類似性**

Belongs to the nucleoplasmin family.


**ドメイン**

The acidic tract A2 mediates histone binding.

**細胞内局在**

Nucleus. Found in the oocyte nucleus before nuclear membrane breakdown, after which it is redistributed to the cytoplasm.

**画像**



**All lanes :** Anti-NPM2 antibody (ab67403) at 1/500 dilution

**All lanes :** NPM2 transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Secondary**

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size:** 24 kDa

**Observed band size:** 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors