

Product datasheet

Anti-ASCC1 antibody ab67919

画像数 1

製品の概要

製品名	Anti-ASCC1 antibody
製品の詳細	Mouse polyclonal to ASCC1
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	ASCC1 (NP_057031.2, 1 a.a. ~ 357 a.a) full-length human protein.
ポジティブ・コントロール	ASCC1 transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67919** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 46 kDa (predicted molecular weight: 46 kDa).

ターゲット情報

機能	Enhances NF-kappa-B, SRF and AP1 transactivation.
組織特異性	Ubiquitous.
配列類似性	Contains 1 KH domain.

## Anti-ASCC1 antibody 画像



Western blot - ASCC1 antibody (ab67919)

**All lanes** : Anti-ASCC1 antibody (ab67919)  
at 1/500 dilution

**Lane 1** : ASCC1 transfected 293T cell lysate  
at 25 µg

**Lane 2** : Non transfected 293T cell lysate at  
25 µg

### Secondary

Goat anti-mouse IgG (H&L)-HRP conjugate at  
1/2500 dilution

**Predicted band size** : 46 kDa

**Observed band size** : 46 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors