

## Product datasheet

# Sheep Anti-Rabbit IgG H&L ab6797

### 製品の概要

製品名	Sheep Anti-Rabbit IgG H&L
由来種	Sheep
ターゲット生物種	Rabbit
アプリケーション	<b>適用あり:</b> IP, Immunodiffusion, Conjugation, ICC/IF, IHC-P, IHC-Fr, WB, ELISA
免疫原	Rabbit IgG whole molecule

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 0.424% Potassium phosphate
精製度	IgG fraction
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads
ポリ/モノ	ポリクローナル

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab6797** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.
ICC/IF		1/2000 - 1/10000.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
WB		1/4000 - 1/20000.
ELISA		1/40000 - 1/200000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors