

## Product datasheet

# Anti-DDX43 antibody ab68029

### 画像数 1

#### 製品の概要

製品名	Anti-DDX43 antibody
製品の詳細	Mouse polyclonal to DDX43
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-648 of Human DDX43
ポジティブ・コントロール	Transfected 293T cell lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.40 Constituent: PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab68029** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

#### 追加情報

WB: 1/500 - 1/1000. Detects a band of approximately 73 kDa (predicted molecular weight: 73 kDa).  
Detection pair

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


### 組織特異性

Expressed in testis. Expressed in many tumors of various histological types at a level that is 100-fold higher than the level observed in normal tissues except testis.

### 配列類似性

Belongs to the DEAD box helicase family.  
Contains 1 helicase ATP-binding domain.  
Contains 1 helicase C-terminal domain.  
Contains 1 KH domain.

## 画像



**All lanes** : Anti-DDX43 antibody (ab68029) at 1/500 dilution

**Lane 1** : DDX43 transfected 293T cell lysate

**Lane 2** : non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size:** 73 kDa

**Observed band size:** 73 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors