

Product datasheet

Anti-PCBD2 antibody ab68113

画像数 1

製品の概要

| | |
|----------|--|
| 製品名 | Anti-PCBD2 antibody |
| 製品の詳細 | Mouse polyclonal to PCBD2 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Human |
| 免疫原 | PCBD2 (NP_115527.2, 1 a.a. ~ 129 a.a) full-length human protein. |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Preservative: None Constituents: Whole serum |
| 精製度 | Whole antiserum |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab68113** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | 1/500 - 1/1000. Detects a band of approximately 14 kDa (predicted molecular weight: 12 kDa). |

ターゲット情報

機能 Involved in tetrahydrobiopterin biosynthesis. Seems to both prevent the formation of 7-pterins


and accelerate the formation of quinonoid-BH2.

Regulates the dimerization of homeodomain protein HNF-1-alpha and enhances its transcriptional activity.

配列類似性

Belongs to the pterin-4-alpha-carbinolamine dehydratase family.

画像



Western blot - PCBD2 antibody (ab68113)

All lanes : Anti-PCBD2 antibody (ab68113)

at 1/500 dilution

Lane 1 : PCBD2 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP

Conjugate secondary

antibody at 1/2500 dilution

Developed using the ECL technique.

Predicted band size: 12 kDa

Observed band size: 14 kDa

Additional bands at: 12 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors