

## Product datasheet

# Anti-KLHDC4 antibody ab68210

### 画像数 1

#### 製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-KLHDC4 antibody  |
| 製品の詳細        | Mouse polyclonal to KLHDC4  |
| アプリケーション     | 適用あり: WB  |
| 種交差性         | 交差種: Human  |
| 免疫原          | Full length protein, corresponding to amino acids 1-520 of Human KLHDC4 |
| ポジティブ・コントロール | KLHDC4 transfected 293T cell lysate                                     |

#### 製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Preservative: None<br>Constituents: PBS, pH 7.2   |
| 精製度    | Protein A purified  |
| ポリ/モノ  | ポリクローナル   |
| アイソタイプ | IgG   |

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab68210** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| WB       |           | 1/500 - 1/1000. Detects a band of approximately 70 kDa (predicted molecular weight: 58 kDa). |

#### ターゲット情報

|       |                           |
|-------|---------------------------|
| 配列類似性 | Contains 6 Kelch repeats. |
|-------|---------------------------|

## 画像



**All lanes** : Anti-KLHDC4 antibody (ab68210)  
at 1/500 dilution

**Lane 1** : KLHDC4 transfected 293T cell  
lysate

**Lane 2** : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at  
1/2500 dilution

**Predicted band size** : 58 kDa

**Observed band size** : 70 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors