

Product datasheet

Anti-RBM6 antibody ab68297

画像数 1

製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-RBM6 antibody  |
| 製品の詳細        | Mouse polyclonal to RBM6  |
| 由来種          | Mouse   |
| 特異性          | This antibody is reactive against mammalian transfected lysate. |
| アプリケーション     | <b>適用あり:</b> WB   |
| 種交差性         | <b>交差種:</b> Human   |
| 免疫原          | RBM6 (AAH46643.1, 1 a.a. ~ 601 a.a) full-length human protein.  |
| ポジティブ・コントロール | RBM6 transfected 293T cell lysate                               |

製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Preservative: None<br>Constituents: 1X PBS, pH 7.2  |
| 精製度    | Protein A purified  |
| ポリ/モノ  | ポリクローナル   |
| アイソタイプ | IgG   |

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab68297** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項  |
|----------|-----------|---|
| WB       |           | 1/500 - 1/1000. Detects a band of approximately 70 kDa. |

|       |  |
|-------|--|
| 機能    | Specifically binds poly(G) RNA homopolymers in vitro.                        |
| 組織特異性 | Ubiquitous in adults.  |
| 配列類似性 | Contains 1 G-patch domain.<br>Contains 1 RRM (RNA recognition motif) domain. |
| 細胞内局在 | Nucleus.   |

## 画像



Western blot - RBM6 antibody (ab68297)

**All lanes :** Anti-RBM6 antibody (ab68297) at 1/500 dilution

**Lane 1 :** RBM6 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Observed band size:** 70 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors