

## Product datasheet

# Anti-AP3M2 antibody ab68535

### 画像数 1

#### 製品の概要

製品名	Anti-AP3M2 antibody
製品の詳細	Mouse polyclonal to AP3M2
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Full length protein (Human, aa 1 to ~418)
ポジティブ・コントロール	AP3M2-transfected 293T cell lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab68535** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 40 kDa (predicted molecular weight: 47 kDa).

#### ターゲット情報


**関連性** Part of the AP-3 complex, AP3M3 is an adaptor-related complex not clathrin-associated. The

complex is associated with the Golgi region and to more peripheral structures. It facilitates the budding of vesicles from the Golgi membrane and may be directly involved in trafficking to lysosomes.

細胞内局在

Golgi apparatus.

画像



Western blot - AP3M2 antibody (ab68535)

**All lanes** : Anti-AP3M2 antibody (ab68535) at 1/500 dilution

**Lane 1** : AP3M2-transfected 293T cell lysate

**Lane 2** : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

**Predicted band size** : 47 kDa

**Observed band size** : 40 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors