

Product datasheet

Anti-REPS1 antibody ab69221

画像数 1

製品の概要

製品名	Anti-REPS1 antibody
製品の詳細	Mouse polyclonal to REPS1
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human protein (entry AAH21211, aa 1-653)
ポジティブ・コントロール	REPS1-transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab69221** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 90 kDa (predicted molecular weight: 81 kDa).

ターゲット情報

機能	May coordinate the cellular actions of activated EGF receptors and Ral-GTPases.
配列類似性	Contains 1 EF-hand domain. Contains 2 EH domains.
翻訳後修飾	EGF stimulates phosphorylation on Tyr-residues.
細胞内局在	Membrane > clathrin-coated pit. Colocalize with ITSN1 at the plasma membrane in structures that are most probably clathrin-coated pits.

画像



All lanes : Anti-REPS1 antibody (ab69221) at 1/500 dilution

Lane 1 : REPS1 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

Predicted band size: 81 kDa

Observed band size: 90 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors