

## Product datasheet

# Anti-NFYC antibody ab69625

### 画像数 1

#### 製品の概要

製品名	Anti-NFYC antibody
製品の詳細	Rabbit polyclonal to NFYC
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat
免疫原	Synthesized peptide derived from the N terminus of human NFYC.
ポジティブ・コントロール	<div style="border: 1px solid #ccc; padding: 5px; display: inline-block;"> <a href="#">ウェスタン・ブロット用ポジティブコントロール Recombinant Human NFYC protein</a> </div> Extracts from 293 cells and HeLa cells

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab69625** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 50 kDa (predicted molecular weight: 50 kDa).

アプリケーション	Abreviews	特記事項
ELISA		1/10000.

## ターゲット情報

機能	Stimulates the transcription of various genes by recognizing and binding to a CCAAT motif in promoters, for example in type 1 collagen, albumin and beta-actin genes.
配列類似性	Belongs to the NFYC/HAP5 subunit family.
細胞内局在	Nucleus.

## 画像



**All lanes :** Anti-NFYC antibody (ab69625) at 1/500 dilution

**Lane 1 :** Extracts from 293 cells with no immunizing peptide

**Lane 2 :** Extracts from HeLa cells with no immunizing peptide

**Lane 3 :** Extracts from HeLa cells with immunizing peptide at 5 µg

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 50 kDa

**Observed band size:** 50 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors