

Product datasheet

Anti-ZFAND1 antibody ab69664

画像数 1

製品の概要

製品名	Anti-ZFAND1 antibody
製品の詳細	Mouse polyclonal to ZFAND1
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human ZFAND1 protein (AAH56877)
ポジティブ・コントロール	ZFAND1 transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab69664** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

追加情報

WB: 1/500 - 1/1000. Detects a band of approximately 18 kDa (predicted molecular weight: Isoform 2, 18 kDa; Isoform 1, 31 kDa).

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

配列類似性

Contains 2 AN1-type zinc fingers.

画像



Western blot - ZFAND1 antibody (ab69664)

All lanes : Anti-ZFAND1 antibody (ab69664)
at 1/500 dilution

Lane 1 : ZFAND1 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at
1/2500 dilution

Predicted band size : 18 kDa

Observed band size : 18 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors