

Product datasheet

Anti-TROAP antibody ab70105

画像数 1

製品の概要

製品名	Anti-TROAP antibody
製品の詳細	Mouse polyclonal to TROAP
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	TROAP (AAH32583.1, 1 a.a. ~ 778 a.a) full-length human protein.
ポジティブ・コントロール	TROAP transfected 293T lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab70105** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 100 kDa (predicted molecular weight: 84 kDa).

ターゲット情報

機能	Could be involved with bystin and trophinin in a cell adhesion molecule complex that mediates an initial attachment of the blastocyst to uterine epithelial cells at the time of the embryo implantation.
組織特異性	Strong expression at implantation sites. Was exclusively localized to the apical side of the syncytiotrophoblast. Also found in macrophages.
細胞内局在	Cytoplasm.

画像



All lanes : Anti-TROAP antibody (ab70105)
at 1/500 dilution

Lane 1 : TROAP transfected 293T lysate

Lane 2 : Non-transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP
at 1/2500 dilution

Predicted band size: 84 kDa

Observed band size: 100 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors