

Product datasheet

Anti-STK38 antibody ab71216

画像数 1

製品の概要

製品名	Anti-STK38 antibody
製品の詳細	Rabbit polyclonal to STK38
由来種	Rabbit
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human
免疫原	Synthetic peptide (Human) from the central region of STK38, conjugated to KLH
ポジティブ・コントロール	mouse brain tissue lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab71216** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/1000.
WB		1/100 - 1/500. Detects a band of approximately 54 kDa (predicted molecular weight: 54 kDa).

ターゲット情報

機能	Negative regulator of MAP3K1/2 signaling. Converts MAP3K2 from its phosphorylated form to its nonphosphorylated form and inhibits autophosphorylation of MAP3K2.
組織特異性	Ubiquitously expressed with highest levels observed in peripheral blood leukocytes.
配列類似性	Belongs to the protein kinase superfamily, AGC Ser/Thr protein kinase family. Contains 1 AGC-kinase C-terminal domain. Contains 1 protein kinase domain.
翻訳後修飾	ISGylated.
細胞内局在	Nucleus. Cytoplasm.

## 画像



Anti-STK38 antibody (ab71216) at 1/60 dilution + mouse brain tissue lysate at 35 µg

**Predicted band size:** 54 kDa

**Observed band size:** 54 kDa

**Additional bands at:** 30 kDa, 40 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors