

## Product datasheet

# Anti-RHBDL2 antibody ab71866

### 画像数 1

#### 製品の概要

製品名	Anti-RHBDL2 antibody
製品の詳細	Mouse polyclonal to RHBDL2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-121 of Human RHBDL2 (AAH13103.1, or Swissprot id Q9NX52-3 without the first 12 amino acids)
ポジティブ・コントロール	RHBDL2 transfected 293T lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab71866** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 12 kDa (predicted molecular weight: 14 kDa).

#### ターゲット情報

**関連性** RHBDL2, rhomboid like protein 2, is involved in regulated intramembrane proteolysis and the subsequent release of functional polypeptides from their membrane anchors. There are 2 named isoforms produced by alternative splicing.

**細胞内局在** Membrane; Multi pass membrane protein.

## 画像



**All lanes** : Anti-RHBDL2 antibody (ab71866)  
at 1/500 dilution

**Lane 1** : RHBDL2 transfected 293T lysate

**Lane 2** : Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP  
at 1/2500 dilution

**Predicted band size:** 14 kDa

**Observed band size:** 12 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors