


Product datasheet

Anti-FAIM2 antibody ab72113

画像数 2

製品の概要

製品名	Anti-FAIM2 antibody
製品の詳細	Rabbit polyclonal to FAIM2
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Mouse, Rat, Dog, Human, Monkey 交差が予測される動物種: Chicken, Opossum 
免疫原	A portion of amino acids 1-50 of human FAIM2 isoform A.
ポジティブ・コントロール	Human, mouse and rat brain lysates. Human cerebellum tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium azide Constituents: 0.2% Gelatin, PBS
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72113** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 3 µg/ml. Detects a band of approximately 29 kDa (predicted molecular weight: 35 kDa).

アプリケーション	Abreviews	特記事項
----------	-----------	------


IHC-P

Use a concentration of 10 µg/ml.

## ターゲット情報

機能	Protects cells uniquely from Fas-induced apoptosis.
組織特異性	Widely expressed. Expressed at high levels in the brain especially in the hippocampus.
配列類似性	Belongs to the B11 family.
細胞内局在	Cell membrane. Cell junction > synapse > postsynaptic cell membrane.

## 画像



**Lanes 1-2 :** Anti-FAIM2 antibody (ab72113)

at 2 µg/ml

**Lane 3 :** Anti-FAIM2 antibody (ab72113) at 1

µg/ml

**Lane 4 :** Anti-FAIM2 antibody (ab72113) at

0.5 µg/ml


**Lane 1 :** Human brain tissue lysate

**Lane 2 :** Human brain tissue lysate with  
Immunizing peptide

**Lane 3 :** Mouse brain tissue lysate

**Lane 4 :** Rat brain tissue lysate

**Predicted band size:** 35 kDa



Analysis of FAIM2 expression in human

cerebellum tissue, using ab72113 at 10 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors