

Product datasheet

Anti-ZFP37 antibody ab72714

画像数 2

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-ZFP37 antibody |
| 製品の詳細 | Mouse polyclonal to ZFP37 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB, ICC/IF |
| 種交差性 | 交差種: Human |
| 免疫原 | Full length human ZFP37 protein (NP_003399.1, 1 a.a. ~ 630 a.a) |
| ポジティブ・コントロール | ZFP37 transfected 293T cell lysate. HeLa cells. |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | pH: 7.20 Constituent: PBS |
| 精製度 | Protein A purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72714** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| WB | | |
| ICC/IF | | |

追加情報

ICC/IF: Use at a concentration of 10 µg/ml. Note: The antibody was purified with Protein A prior to use in immunofluorescence.

WB: 1/500 - 1/1000. Detects a band of approximately 80 kDa (predicted molecular weight: 68 kDa).

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性 ZFP37 may be involved in transcriptional regulation.

細胞内局在 Nuclear

画像



Western blot - ZFP37 antibody (ab72714)

All lanes : Anti-ZFP37 antibody (ab72714) at 1/500 dilution

Lane 1 : ZFP37 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate


Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 68 kDa

Observed band size: 80 kDa



Immunocytochemistry/ Immunofluorescence - ZFP37 antibody (ab72714)

Immunofluorescence analysis of ZFP37 expression in HeLa cells using purified ab72714 at 10 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors