

Product datasheet

Anti-BAT4 antibody ab72716

1 References [画像数 1](#)

製品の概要

製品名	Anti-BAT4 antibody
製品の詳細	Mouse polyclonal to BAT4
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human BAT4 protein (NP_149417, 1 a.a. ~ 356 a.a)
ポジティブ・コントロール	BAT4 transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72716** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 35 kDa (predicted molecular weight: 39 kDa).

ターゲット情報


関連性 BAT4 (HLA-B-associated transcript 4) contains two ANK repeats and one G-patch domain. It is

thought to be involved in some aspects of immunity.

細胞内局在

Endoplasmic reticulum

画像



Western blot - BAT4 antibody (ab72716)

All lanes : Anti-BAT4 antibody (ab72716) at 1/500 dilution

Lane 1 : BAT4 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size : 39 kDa

Observed band size : 35 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors