

## Product datasheet

# Anti-CPA5 antibody ab72935

### 画像数 1

#### 製品の概要

製品名	Anti-CPA5 antibody
製品の詳細	Mouse polyclonal to CPA5
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human CPA5 protein (NP_525124)
ポジティブ・コントロール	CPA5 transfected 293T cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72935** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 49 kDa.

#### ターゲット情報

組織特異性

Expression is very low or not detectable.


配列類似性

Belongs to the peptidase M14 family.

細胞内局在

Secreted.

## 画像



**All lanes** : Anti-CPA5 antibody (ab72935) at 1/500 dilution

**Lane 1** : CPA5 transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/5000 dilution

Developed using the ECL technique.

**Predicted band size:** 49 kDa

**Observed band size:** 49 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors