

## Product datasheet

### Anti-C19orf46 antibody ab73050

#### 画像数 1

#### 製品の概要

製品名	Anti-C19orf46 antibody
製品の詳細	Mouse polyclonal to C19orf46
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-291 of Human C19orf46
ポジティブ・コントロール	Transfected 293T cell line lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: Whole serum
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab73050** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

#### 追加情報

WB: 1/500 - 1/1000. Expected band of 32 kDa (Isoform 2). Detects a band of approximately 37 kDa.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

### 細胞内局在

Nucleus membrane; Single-pass type IV membrane protein

## 画像



**All lanes** : Anti-C19orf46 antibody (ab73050)  
at 1/500 dilution

**Lane 1** : C19orf46 transfected 293T cell  
lysate

**Lane 2** : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat anti-mouse IgG (H&L)-HRP  
conjugate at 1/2500 dilution

**Observed band size:** 37 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors