


Product datasheet

Anti-DDX55 antibody ab74816

画像数 1

製品の概要

製品名	Anti-DDX55 antibody
製品の詳細	Rabbit polyclonal to DDX55
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Mouse 
免疫原	Synthetic peptide derived from an internal sequence of human DDX55
ポジティブ・コントロール	extracts from HeLa cells.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab74816** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 68 kDa (predicted molecular weight: 68 kDa).
ELISA		1/40000.

ターゲット情報

機能	Probable ATP-binding RNA helicase.
配列類似性	Belongs to the DEAD box helicase family. DDX55/SPB4 subfamily. Contains 1 helicase ATP-binding domain. Contains 1 helicase C-terminal domain.
ドメイン	The Q motif is unique to and characteristic of the DEAD box family of RNA helicases and controls ATP binding and hydrolysis.

画像



All lanes : Anti-DDX55 antibody (ab74816) at 1/500 dilution

Lane 1 : extracts from HeLa cells

Lane 2 : extracts from HeLa cells with immunising peptide at 5 µg

Lysates/proteins at 5 µg per lane.

Predicted band size: 68 kDa

Observed band size: 68 kDa

Additional bands at: 31 kDa, 80 kDa, 85 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors