

Product datasheet

Anti-ARF4 antibody ab75073

★☆☆☆☆ 2 Abreviews [画像数 1](#)

製品の概要

製品名	Anti-ARF4 antibody
製品の詳細	Rabbit polyclonal to ARF4
由来種	Rabbit
特異性	ARF4 antibody detects endogenous levels of total ARF4 protein.
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Synthetic peptide (Human)
ポジティブ・コントロール	LOVO cell lysate, K562 cell lysate, HT-29 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
精製度	Immunogen affinity purified
特記事項 (精製)	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab75073** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★☆☆☆☆	1/500 - 1/1000. Predicted molecular weight: 21 kDa.
ELISA		1/40000.

ターゲット情報

機能	GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.
配列類似性	Belongs to the small GTPase superfamily. Arf family.
細胞内局在	Golgi apparatus.

画像



All lanes : Anti-ARF4 antibody (ab75073) at 1/1000 dilution

Lane 1 : LOVO cell lysate

Lane 2 : K562 cell lysate

Lane 3 : HT-29 cell lysate

Lane 4 : K562 cell lysate with Blocking peptide at 5 µg

Lysates/proteins at 5 µg per lane.

Predicted band size: 21 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors