

Product datasheet

Anti-zinc finger protein 138 antibody ab76445

画像数 1

製品の概要

製品名	Anti-zinc finger protein 138 antibody
製品の詳細	Mouse monoclonal to zinc finger protein 138
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Recombinant fragment 交差が予測される動物種: Human
免疫原	Recombinant fragment: TNLSKPKKIH TGEKPYKCEV CGKAFHQSSI LTKHKIIRTG EKPYPKCAHCG KAFKQSSHLT RHKIIHTEEK PYKCEQCGKV FKQSPTLTKH QIYTGEEPYK, corresponding to amino acids 151-251 of Human zinc finger protein 138, with a 26kDa enzyme tag. Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76445** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

ELISA

追加情報

ELISA: Use at an assay dependent dilution.

WB: 1/500 - 1/1000. Detects a band of approximately 40 kDa.

Note: This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報


関連性

May be involved in transcriptional regulation as a repressor.

細胞内局在

Nuclear

画像



Anti-zinc finger protein 138 antibody
(ab76445) at 1/500 dilution + immunogen at
0.2 µg

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at
1/5000 dilution

Observed band size : 40 kDa

Western blot - zinc finger protein 138 antibody
[4D11] (ab76445)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors