

## Product datasheet

# Anti-REEP5 antibody ab76451

★★★★☆ 2 Abreviews [画像数 1](#)

### 製品の概要

製品名	Anti-REEP5 antibody
製品の詳細	Mouse monoclonal to REEP5
アプリケーション	<b>適用あり:</b> WB, ELISA, ICC/IF
種交差性	<b>交差種:</b> Human, Xenopus laevis
免疫原	Recombinant fragment: KCGFLLWCMA PSPSNGAELL YKRIIRPFFL KHESQMDSVV KDLKDKAKET ADAITKEAKK ATVNLLGEEK KS, corresponding to amino acids 113-185 of Human REEP5, with tag <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76451** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★★★★☆	
ELISA		

ICC/IF



## 追加情報

ELISA: Use at an assay dependent dilution.

WB: 1/500 - 1/1000. Predicted molecular weight: 21 kDa.

Note: This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


## 関連性

May promote functional cell surface expression of olfactory receptors.

## 細胞内局在

Cell Membrane

## 画像



Western blot - REEP5 antibody [3G11] (ab76451)

Anti-REEP5 antibody (ab76451) at 1/500 dilution + Recombinant protein at 0.2 µg

**Secondary**

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/5000 dilution

**Predicted band size** : 21 kDa

**Observed band size** : 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors