

Product datasheet

# Anti-Polynuclear Aromatic Hydrocarbons antibody [BAP-13] ab768

4 References

製品の概要

製品名	Anti-Polynuclear Aromatic Hydrocarbons antibody [BAP-13]
製品の詳細	Mouse monoclonal [BAP-13] to Polynuclear Aromatic Hydrocarbons
由来種	Mouse
特異性	Tested compound : Reactivity ± SD [%]: Benzo[a]pyrene = 100 Naphthalene = <1 Acenaphthene = 1 ± 0.2 Acenaphthylene = <1 Fluorene = 5 ± 1 Phenanthrene = 23 ± 4 Anthracene = 25 ± 4 Fluoranthene = 60 ± 9 Pyrene = 47 ± 7 Benzo[a]anthracene = 44 ± 7 Chrysene = 42 ± 6 Benzo[b]fluoranthene = 26 ± 4 Benzo[k]fluoranthene = 15 ± 2 Dibenzo[a,h]anthracene = 9 ± 1 Benzo[g,h,i]perylene = 15 ± 2 Indeno[1,2,3-c,d]pyrene = 63 ± 10 Benzo[e]pyrene = 55 ± 9 Perylene = 8 ± 1
アプリケーション	<b>適用あり:</b> ELISA
免疫原	Benzo[a]pyrenyl-1-butyric acid conjugated to BSA.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 15mM Sodium Azide Constituents: PBS, pH 7.4
精製度	>95% by SDS-PAGE
特記事項 (精製)	Purified from ascites using precipitation methods and DEAE-chromatography.
ポリ/モノ	モノクローナル
クローン名	BAP-13
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab768** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

ELISA

#### 追加情報

ELISA: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

#### ターゲット情報

#### 関連性

Polynuclear aromatic hydrocarbons or Polycyclic Aromatic Hydrocarbons (PAHs) are contain multiple benzene rings. They are generally found as components of fuels, oils and greases.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors