

Product datasheet

Anti-PPCDC antibody ab76903

画像数 1

製品の概要

製品名	Anti-PPCDC antibody
製品の詳細	Mouse polyclonal to PPCDC
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human PPCDC protein (AAH14409.1).
ポジティブ・コントロール	PPCDC transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76903** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 22 kDa.

ターゲット情報

機能	Necessary for the biosynthesis of coenzyme A. Catalyzes the decarboxylation of 4-phosphopantothenoylcysteine to form 4'-phosphopantotheine.
パスウェイ	Cofactor biosynthesis; coenzyme A biosynthesis; CoA from (R)-pantothenate: step 3/5.

## Anti-PPCDC antibody 画像



**All lanes** : Anti-PPCDC antibody (ab76903)  
at 5 µg/ml

**Lane 1** : PPCDC transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Secondary**

Goat Anti-Mouse IgG (H&L)-HRP at 1/2500  
dilution

**Predicted band size** : 22 kDa

**Observed band size** : 22 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors