

Product datasheet

Anti-GHRH antibody ab80008

画像数 1

製品の概要

製品名	Anti-GHRH antibody
製品の詳細	Rabbit polyclonal to GHRH
由来種	Rabbit
アプリケーション	適用あり: ELISA, IHC-P
種交差性	交差種: Human
免疫原	KLH conjugated synthetic peptide selected from the N terminal region of Human GHRH.
ポジティブ・コントロール	Human lung carcinoma tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Ammonium Sulphate Precipitation
特記事項 (精製)	This antibody is prepared by Saturated Ammonium Sulfate precipitation followed by dialysis against PBS.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab80008** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/1000.
IHC-P		1/10 - 1/50.

ターゲット情報

機能	GRF is released by the hypothalamus and acts on the adenohypophyse to stimulate the secretion of growth hormone.
配列類似性	Belongs to the glucagon family.
細胞内局在	Secreted.

画像



ab80008, at a 1/10 dilution, staining GHRH in formalin fixed, paraffin embedded human lung carcinoma tissue by Immunohistochemistry. ab80008 was peroxidase conjugated to the secondary antibody, followed by DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - GHRH antibody - N-terminal (ab80008)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors