

Product datasheet

Anti-Corticotropin Releasing Factor antibody ab80359

画像数 1

製品の概要

製品名	Anti-Corticotropin Releasing Factor antibody
製品の詳細	Chicken polyclonal to Corticotropin Releasing Factor
由来種	Chicken
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to amino acids 13-87 of human Corticotropin Releasing Factor derived from <i>E. coli</i>
ポジティブ・コントロール	<i>E. coli</i> derived Corticotropin Releasing Factor fusion protein.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab80359** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/3000. Predicted molecular weight: 21 kDa.

ターゲット情報

関連性 CRF precursor is a 41-amino acid peptide derived from a 196-amino acid preprohormone. CRF is a neuroendocrine hormone secreted by the paraventricular nucleus (PVN) of the hypothalamus and is a key player in the hypothalamic pituitary adrenal axis. It is also synthesized in peripheral tissues, such as T lymphocytes and is highly expressed in the placenta. Neuronal CRF is important in a variety of diverse roles including, memory, depression, stress, kidney function, parturition, metabolism and feeding (to mention a few).

細胞内局在 Cytoplasmic and Secreted

画像



Anti-Corticotropin Releasing Factor antibody (ab80359) at 1/5000 dilution + *E. coli* derived Corticotropin Releasing Factor fusion protein at 0.4 µg

Secondary

goat-anti-IgY/Fc-HRP at 1/10000 dilution

Developed using the ECL technique.

Predicted band size: 21 kDa

Observed band size: 21 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors